StatusBlue

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Instructions for use English

Product description

StatusBlue is a diagnostic impression material based on A-silicones. The material is available in the Automix Cartridge and in the ready-to-use pre-filled MixStar Cartridge, which is suitable for use in the MixStar eMotion and other commercially available mixing machines.

Intended purpose

- Impressions for provisional repairs

- Situation Impressions
- Opposing jaw models
- Impressions for the fabrication of vacuum-formed stents

Limitations of use

There are no known limitations of use.

Indication

All dental situations requiring diagnostic casts e.g. tooth loss, caries and trauma.

Contraindication

Do not use the material where there are known allergies to any of the components or contact allergies.

Patient target group

Persons receiving treatment as part of a dental procedure.

Intended user

Dentist, dental assistant staff

Notes for use

- The recommended times can be found in the chapter headed “Timing.”
- In order to achieve an optimum mixing result, DMG recommends the use of the mixing tips available from DMG. All mixing tips are available as refills.
- StatusBlue can be used in all common impression trays.
- In order to improve the adhesion of the material to metal or plastic trays, a tray adhesive can be applied.

Using the MixStar Cartridge in the MixStar eMotion (Fig. see “Handling”)

Follow the instructions for use for the “MixStar eMotion” unit!

1. Slide the locking mechanism upwards in order to open the MixStar Cartridge. Do not remove the locking mechanism from the MixStar Cartridge [Fig. 1].
2. Remove the seal from the MixStar Cartridge and discard. Attach a new mixing tip [Fig. 2].
3. Secure the mixing tip by sliding the locking mechanism into the locked position [Fig. 3].
4. Insert the MixStar Cartridge into the MixStar eMotion [Fig. 4]. A brief acoustic signal confirms correct insertion of the cartridge.

Note: Prior to using a new cartridge for the first time, discard approx. the first 3 cm of the ejected material.

Note: After the application, leave the used mixing tip on the MixStar Cartridge as a seal.

Using the MixStar Cartridge in mixing units from other manufacturers

Follow the instructions for use for the mixing unit used!

- First insert the MixStar Cartridge into the mixing unit without the mixing tip. Otherwise appropriate mixing of the material cannot be guaranteed!
- In this case proceed as described in the section “Changing the mixing tip.”

Changing the mixing tip

Follow the instructions for use for the mixing unit used!

Note: Attach a new mixing tip each time before mixing.

1. Slide the locking mechanism upwards.
2. Remove the old mixing tip as well as any hardened material blocking the cartridge openings.
3. Place the mixing tip securely on the corresponding MixStar Cartridge opening.
4. Lock the mixing tip in place with the locking mechanism.

Using the Automix Cartridge (Fig. see “Handling”)

Using the Automix Cartridge

1. Push up the lever [A] at the back of the Automix Dispenser and pull back the slide [B] all the way.
2. Lift up the plastic lever [C], insert the Automix Cartridge and lock it using the plastic lever [C].

Note: Make sure that the notches on the Automix Cartridge line up with those on the Automix Dispenser.

Attaching a mixing tip

1. Remove the cartridge cap [D] or used mixing tip after turning it 90° counterclockwise, and discard it.
2. Attach a new mixing tip [E].

Note: In order to achieve an optimum mixing result, DMG recommends the use of the mixing tips available from DMG. All mixing tips are also available in refill packs.

Note: Make sure that the notches on the mixing tip line up with those on the Automix Cartridge.

3. Twist the mixing tip clockwise by 90° to lock it.

Applying the material

- When using an Automix Cartridge for the first time, squeeze out a pea-sized amount of material and discard it.
- Mix the material in the mixing tip using the lever [F] on the Automix Dispenser and apply directly.

Note: After the application, leave the used mixing tip on the Automix Cartridge as a seal!

Removing the Automix Cartridge

1. Push up the lever [A] at the back of the Automix Dispenser and pull back the slide [B] all the way.
2. Lift up the plastic lever [C] and remove the Automix Cartridge.

Disinfecting the Automix Dispenser

- The Automix Dispenser can be reused after disinfection, although it should be replaced if there are visible signs of wear.
- Disinfection can be carried out using a commercially available dipping disinfectant. DMG recommends using only disinfectant solutions listed by the RKI (Robert Koch Institute).
- The Automix Dispenser is also autoclavable.

Recommended use

1. Fill the prepared tray with the material. When doing so, hold the impression tray at a slight angle beneath the mixing tip.

Note: To avoid air bubbles, the end of the mixing tip should remain immersed in the material throughout application.

2. Position the filled impression tray slowly and gently in the mouth within the working time.
3. Remove the impression tray from the mouth once the setting time has ended (see “Removal of the impression”).

Removal of the impression

- In the area where the teeth have been prepared, remove the impression parallel to the longitudinal axis of the teeth.
- After removal, proceed as described in the section “Cleaning and disinfection of the final impression.”

Cleaning and disinfection of the final impression

- Only disinfect the final impression using dipping solutions.
- DMG recommends the use of a 2 % glutaraldehyde solution.
- Ensure that the disinfectant used is indicated for use with addition-cured silicone impression material.
- Follow the manufacturer’s instructions for the disinfectant solution precisely!
- Leave the final impression to dry and store until casting (see “Storage of the impression”).
- For more information, see the statement issued by the German Society for Dental and Oral Medicine (DGZMK). Comments of the DGZMK 3/93 V 1.0; www.dgzmk.de
- ADA recommends following the advice of the Council on Dental Materials, Instruments and Equipment. JADA 1991;122(3):110; www.ada.org

Casting and galvanization

- The impression should be cast 10 min after removal from the mouth at the earliest.
- The use of a special purpose hard plaster (ISO 6873:2013, Type 3; e.g. Moldano, Bayer) is recommended for the cast and an ultra hard dental plaster for the models (ISO 6873:2013, Type 4; e.g. Fujirock II, GC)
- The impression can be copper-plated or silver-plated.

Storing the impression

- Store the impression at 15 – 25 °C (59 – 77 °F) in a dry place!
- Keep away from direct sunlight.
- The impression is stable and can be cast even after months without loss of accuracy.

Timing

StatusBlue		
Working time including mixing time * [s]	≤ 75	≤ 75
Setting time in the mouth [s]	≥ 105	≥ 105

StatusBlue		
------------	--	--

* The specified times are valid at a room temperature of 23 °C and a normal relative humidity of 50 %. Higher temperatures shorten, lower temperatures lengthen these times. The working time describes the period from when mixing begins to the time by which the impression tray must be positioned in the patient's mouth.

Interactions

- Rinse off any residual retraction solution thoroughly with water, otherwise this may lengthen the setting time in some instances.
- Never allow the impression to come into contact with organic solvents or liquids containing solvents as this may result in a swelling of the material and may therefore lead to inaccuracies in the cast.
- Do not allow the impression material to come into contact with latex gloves. DMG recommends the use of PVC or PE-based gloves.
- Do not mix StatusBlue with condensation curing silicone or polyether impression materials.
- Leftover methacrylate-based materials may alter the setting behavior of the silicone impression materials. If a restoration was prepared beforehand, carefully clean the preparation and the surrounding tissue with alcohol-soaked gauze and water spray before taking the impression.

Residual risks/side effects

- We place great importance on the quality of our products. However, please check the product for damage and visible contamination before use.
- There are no known side effects to date.
- No cases of hypersensitivity or allergic reactions are known to date. However, reactions to components of the material cannot be ruled out.

Warnings/precautions

- For dental use only!
- Avoid contact with the eyes! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.
- Follow the manufacturer's instructions for other products that are used with the material/materials.
- Serious incidents involving this product must be reported to the manufacturer (productsafety@mtcompanies.com) and to the responsible registration authorities.

Product performance characteristics

Meets the requirements of ISO 4823:2025, Type 2: medium-bodied consistency

Storage/disposal

- Store in a dry place at 15 to 25 °C (59 to 77 °F).
- Protect from exposure to direct sunlight!
- After the application, leave the used mixing tip on the cartridge as a seal!
- Do not use after the expiration date!
- Disposal must comply with national regulations.

Composition

StatusBlue (MixStar): SiO₂, polyvinylsiloxane, polydimethylsiloxane, aluminum silicate, methyl hydrogen polysiloxane, additives, pigments, catalysts.

StatusBlue (Automix): SiO₂, polyvinylsiloxane, polydimethylsiloxane, aluminum silicate, methyl hydrogen polysiloxane, additives, pigments, catalysts.

.

Packaging

1 MixStar Cartridge @ 380 ml, 10 MixStar Tips	
	REF 909693

5 MixStar Cartridges @ 380 ml, 50 MixStar Tips	
	REF 909694

2 Automix Cartridges @ 50 ml, 4 Automix Tips	
	REF 909790

8 Automix Cartridges @ 50 ml, 16 Automix Tips	
	REF 909791

Accessories

1 MixStar eMotion 230 V Int.	REF 900480
50 MixStar Tips White	REF 909567
1 Automix Dispenser Type 50 10:1	REF 110411
25 Automix Tips Pink	REF 909797